

OptiSpec® Gynecology Light

Hands-free ultra-bright illumination provides visualization without inconvenience.



Utah Medical Products' patented* OptiSpec Light is a new concept in non-colposcopic illumination of the cervix. An ultra-bright LED selected to provide a pure white light spectrum has been mounted in a small, clip-on disposable package. The result is excellent illumination of the cervix with a device that otherwise seems like it's not even there!

Cervical visualization through a colposcope with a bright white light provides critical visual information with minimal clinician fatigue. However, the use of the colposcope for other every-day exams is impractical. Other methods of cervical illumination emit a dull yellow, low intensity light, and usually require one hand to actively hold the lighting device.



- Compact light clips on to most common vaginal specula
- Unobtrusive configuration remains out of visual and working field
- Provides a simple, hands-free light that improves visualization during:
 - gyn exams
 - pap smears
 - LETZ® procedures
 - diagnosis of abnormal obstetric bleeding
 - ER exams for vaginal trauma
 - any other directly visualized vaginal procedures
- Efficient light-emitting diode (LED) provides truer color visualization with a light that is whiter than halogen bulbs

Order Information

- Item Number: LITE-WS
- Packaging: 25 per box
Single Use, Sterile

Specifications

- Housing: ABS
- Light source: Ultrabright white LED

The OptiSpec Light is LATEX FREE
*U.S. Patent 7,631,981



United States

Utah Medical Products, Inc.
7043 South 300 West
Midvale, Utah 84047
800.533.4984 • 801.566.1200
fax 801.566.2062

Ireland

Utah Medical Products Ltd.
Athlone Business & Technology Park
Athlone, County Westmeath
(090) 647.3932 • fax (090) 647.5608

United Kingdom

Femcare Ltd.
32 Premier Way
Romsey, Hampshire
SO51 9DQ
(01794) 525.100 • fax (01794) 525.101

www.utahmed.com