

Is your first-line technology for detecting endometrial abnormality effective?



Disposable sampling devices obtain a sample from only 4% of the endometrial surface area,¹ and fail to detect 1 in 6 endometrial cancers.²

In-office endometrial sampling is a cost-efficient method for first-line diagnosis of abnormal uterine bleeding (AUB). However, published studies demonstrate that existing suction curette devices do not provide consistent specimen volume or quality, and insinuate the cause may be due to sampling tissue through a single, small port. Consequently, false negative assessment often occurs in patients with focal pathology.

EndoCurette®

Excellent Results for Diagnosis of AUB



EndoCurette uses four bowed curetting elements to remove endometrium independent from the orientation of the four elongated sampling ports. This special configuration is most effective with a single fundus-to-os draw with a twisting motion. EndoCurette obtains robust tissue samples with intact glands and stroma, often providing first-pass definitive diagnosis of abnormal bleeding.

Experience using EndoCurette and see what you've been missing!

1 Rodriguez GC et al, *A comparison of the Pipelle device and the Vabra aspirator as measured by endometrial denudation in hysterectomy specimens: The Pipelle device samples significantly less of the endometrial surface than the Vabra aspirator*, Am J Obstet Gynecol 1993;168:55-9
2 Guido RS et al, *Pipelle endometrial sampling: sensitivity in the detection of endometrial cancer*, J Repr Med 1995;40:553-5

EndoCurette Clear

CUR-120
20/box

EndoCurette Clear's 30 ml syringe maintains vacuum for sampling, even when aspirating fluid.

- Special multi-port tip configuration is designed to obtain a sample representative of a majority of the endometrial surface to improve detection of focal pathology.
- Deliberate fundus-to-os sampling motion prevents patient discomfort and risk of trauma by eliminating repeated fundal contact of tip.
- Tip profile, vacuum plunger, and cannula rigidity provide stiffness that facilitates insertion and may help provide access through mildly stenotic cervix.
- 3mm diameter minimizes patient discomfort.

EndoCurette

CUR-100
25/box

Traditional plunger-cannulus system provides positive vacuum and rigidity during insertion.



United States

Utah Medical Products, Inc.
7043 South 300 West
Midvale, Utah 84047
800.533.4984 • 801.566.1200
fax 801.566.2062
customer@utahmed.com

Ireland

Utah Medical Products Ltd.
Athlone Business & Technology Park
Athlone, County Westmeath N37 XK74
(090) 647.3932 • fax (090) 647.5608
customerserviceireland@utahmed.com

United Kingdom



Femcare Ltd.
32 Premier Way
Romsey, Hampshire SO51 9DQ
(01794) 525.100 • fax (01794) 525.101
enquiries@femcare.co.uk

www.utahmed.com